



# $\alpha$ GEMINI

**NEW**

SCREENING TOOL FOR GLOBAL AROMA/VOC ANALYSIS

## *Simplifying Your Quality Control*

Food

Beverages

Flavors

Fragrances



Perfumes

Cosmetics

Polymers

Chemicals

### *HITTING YOUR TARGET MARKET CONSISTENTLY by*

- Rapidly and objectively controlling your product quality by aroma/VOC fingerprint analysis.
- Verifying incoming raw materials to your gold standard.
- Detecting taints and off-odors before product shipment.
- Push-button system with “ Pass/Fail ” reporting for plant operator level.
- Statistical Quality Control chart to monitor process or quality variability.



MULTI ORGANOLEPTIC SYSTEMS

**THE SPECIALIST IN SMART SENSING SYSTEMS**




# αGEMINI

## GET A Rapid AND OBJECTIVE ANSWER ON PRODUCT OR PROCESS QUALITY !

Out of Specification easily detected

Control chart of product process quality variability



User-defined "Warning Limit" of the product quality variability

Reference and "accepted" unknow product fingerprints

### FEATURES

- ▶ Chemical sensor array technology
- ▶ Sample analysis in less than 5 minutes
- ▶ Almost no sample preparation
- ▶ Supervisor and operator with multi-lingual interface
- ▶ Clear answer : pass or fail
- ▶ Library building ability
- ▶ Data transfer into LIMS
- ▶ High throughput autosampler

### BENEFITS

- ▶ Robustness, low cost maintenance
- ▶ Fast control of aroma/VOC fingerprint
- ▶ Easy implementation
- ▶ Full compliance with factory and quality control laboratory
- ▶ Easy interpretation of results
- ▶ Upgrade the database when needed
- ▶ Tightly supervise quality control
- ▶ Increase plant efficiency and availability

**Alpha M.O.S. France**  
 20, avenue Didier-Daurat  
 31400 Toulouse - France  
 Tél. : +(33) 5.62.47.53.80  
 Fax : +(33) 5.61.54.56.15

**Alpha M.O.S. America**  
 102, Towne Centre Drive  
 Hillsborough - NJ 08876  
 Phone : 908.359.9396  
 Fax : 908.359.9398

Please call me for :

- having some samples analysed
- a quotation
- a discussion on my application
- a presentation of the αGEMINI

Name : .....

Job title : .....

Company : .....

Business area : .....

Address : .....

City / State / Zip : .....

Country : .....

E-mail : .....

Phone : ..... Fax : .....



MULTI ORGANOLEPTIC SYSTEMS

Internet : [www.alpha-mos.com](http://www.alpha-mos.com)  
 North America e-mail : [amusa@alpha-mos.com](mailto:amusa@alpha-mos.com)  
 Europe & Export e-mail : [eur@alpha-mos.com](mailto:eur@alpha-mos.com)